



Model Number		OXYGEN FLOW RATE			OXYGEN PRESSURE		DIMENSIONS		WEIGHT		dB(A)	Power kW
		90%	93%	95%	BAR	PSIG	mm	inches	KG	LBS		
MD-5A	LPM	5.0	4.3	3.8	0.8	12.0	680 L x 445 W x 780 H	26.8 L x 17.5 W x 30.7 H	66	145	≤55	1.1
	SCFH	10.6	9.0	8.0								
MD-10A	LPM	10.0	8.5	7.5	0.8	12.0	682 L x 445 W x 780 H	26.8 L x 17.5 W x 30.7 H	66	145	≤55	1.1
	SCFH	21.2	18.0	15.9								
MD-10B	LPM	10.0	8.5	7.5	4.5	65.0	795 L x 450 W x 795 H	31.3 L x 17.7 W x 31.3 H	102	225	≤55	1.1
	SCFH	21.2	18.0	15.9								
MD-20B	LPM	20.0	17.0	15.0	4.5	65.0	797 L x 450 W x 795 H	31.3 L x 17.7 W x 31.3 H	102	225	≤55	1.8
	SCFH	42.4	36.0	31.8								

*all power supply 110V/ single phase/ 60Hz