

MD series small-scale oxygen generator





Storage Friendly

Whole cabinet structure with casters, can be installed at any place you like

All-In-One Design

All Components In One Cabinet

Super Save Space



Reliable oxygen generation in highland region

The generator produces steady oxygen flow and pressure with purity of $93\pm 3\%$, at 5000 meters above sea level



High precision parts



High air intake



High oxygen recovery rate



Minimized pressure loss



Optimized parameter settings



Advanced gas flow system


All-in-One Small-scale Oxygen Generator

The small-scale oxygen generator producing $93\pm 3\%$ oxygen can be used in animal operations. our machines are suitable for veterinary use of all scales.

Small & Customized

Car power inverter available

Real-time oxygen generation





Type	MD-10	MD-20
Oxygen flowrate(L/min)	10	20
Oxygen purity(%)	93±3%	
Oxygen pressure(bar)	4.0±0.5	
Noise level dB(A)	≤ 55	
L*W*H(mm)	680*445*780	795*450*885
Weight(kg)	66	102
Power supply	AC220V 50Hz/AC110V 60Hz	
Power (kW)	1.1	1.8



Type	MD-10 mini
Oxygen flowrate(L/min)	10
Oxygen purity(%)	93±3%
Oxygen pressure(MPa)	0.4±0.05
Noise level dB(A)	≤ 55
L*W*H(mm)	580*390*590
Weight(kg)	53
Power supply	AC220V 50Hz/AC110V 60Hz
Power (kW)	1.1



Type	MD-50
Oxygen flowrate(L/min)	30~40
Oxygen purity(%)	93±3%
Oxygen pressure(MPa)	0.4±0.05
Noise level dB(A)	≤ 60
L*W*H(mm)	820*585*1820
Weight(kg)	215~245
Power supply	AC220V 50Hz
Power (kW)	3.3~4.2

Product link : <https://www.meditecho2.com/small-scale-oxygen-generators.html>